



In-Tube Extraction

In-Tube Extraction Sample Preparation for Gas Chromatography

...*CTC ANALYTICS*

CTC Analytics AG
Industriestrasse 20
CH-4222 Zwingen
Switzerland

Tel: +41 61 7658100
Fax: +41 61 7658199
E-mail: support@ctc.ch
Web: www.ctc.ch

In-Tube Extraction


ITEX, a new Sample Preparation Technique for Gas Chromatography

Bruno Baltensperger

CTC Analytics AG, Industriestrasse 20, CH-4222 Zwingen, Switzerland

Beat Schilling

BGB Analytik AG, Lettenstrasse 97, CH-8134 Adliswil, Switzerland




Range of possible applications

- Environmental
- Drinking Water
- Foods / Flavours
- Consumer Products
- Forensics / Toxicology
- Petrochemicals
- Polymers
- Pharmaceuticals
- Residual Solvents


Rapid and efficient sample enrichment of volatile and semi-volatile compounds in solid, liquid and gaseous samples

- In-tube extraction and direct thermal desorption using proven industry standard adsorbents
- Syringe only concept for transparent sample handling, no sample loops, transfer lines, or switching valves
- No GC injector modifications, no cryo-focussing necessary
Top mounted on GC's, saves valuable bench space
- Interfaces with any CombiPAL System controlled by all major GC/GC-MS Systems

ITEX Sample Extraction Procedure



Enrichment | Relationship between Sample Volume & Sensitivity



Enrichment of Methyl esters on a Tenax TA ITEX trap. 1 μ l of a mixture of C₄, C₆, C₈, C₁₀ Methyl esters in Methanol (100ng/ μ l) was injected into a 20ml HS-vial. After conditioning at 40°C for 10min. 1ml of the headspace using an ITEX trap without packing material was injected to determine a "static headspace value" (0x). Afterwards the needle was replaced by the TENAX TA ITEX trap. The enrichment of the solutes on the trap was studied using various numbers of pumping strokes.

ITEX Parameter:

Extraction Speed: 100 μ l/sec.

Extraction Strokes: 0 / 10 / 20 / 40 / 60

Temperature Pumping Syringe / Sample Incubation: 40°C/10min.

Desorption at 250°C, 15sec. splitless

Chromatography:


Injection: Splitless 15sec. at 250°C, Carrier gas: 0.2bar Hydrogen

Column: BGB-1 15m x 0.32mm ID, 1.0 μ m film

Temperature Program: 40°C - 1min. - 10°C/min to 200°C

Detection: FID 250°C

Environmental Application | EPA Method 502.2



Comparison of ITEX analysis versus Static Headspace
 Sample: Purge and Trap calibration mix
 (Restek Cat.No. 30431 502.2 CAL2000 Mega-Mix)

Static Headspace Parameter

60°C / 10min / 1ml sample volume

ITEX Parameter:

Extraction Speed: 100µl / sec.

Total Pumping Strokes: 50

Temperature Pumping Syringe / Sample Incubation: 60°C -10min.

Desorption at 200°C, 15sec. splitless

Chromatography:

Injection: Splitless 15sec. at 250°C / Carrier gas: 0.2bar hydrogen


Column: Rtx-502.2 60m x 0.32mm ID, 1.8µm film

Temperature Program: 40°C - 1min. - 10°C - min to 220°C

Detection: FID 250°C

1 1,1-Dichloroethylene	14 Trichloroethylene	28 Ethylbenzene	42 tert-Butylbenzene
2 Methylene chloride (dichloromethane)	15 1,2-Dichloropropane	29 m-Xylene	43 1,2,4-Trimethylbenzene
3 trans 1,2-Dichloroethylene	16 Bromodichloromethane	30 p-Xylene	44 sec-Butylbenzene
4 1,1-Dichloroethane	17 Dibromomethane	31 o-Xylene	45 4-Isopropyloluene (p-Cymene)
5 2,2-Dichloropropane	18 cis-1,3-Dichloropropylene	32 Styrene	46 1,3-Dichlorobenzene
6 cis-1,2-Dichloroethylene	19 Toluene	33 Isopropylbenzene	47 1,4-Dichlorobenzene
7 Chloroform	20 trans-1,3-Dichloropropylene	34 Bromoform	48 n-Butylbenzene
8 Bromochloromethane	21 1,1,2-Trichloroethane	35 1,1,2,2-Tetrachloroethane	49 1,2-Dichlorobenzene
9 1,1,1-Trichloroethane	22 1,3-Dichloropropane	36 1,2,3-Trichloropropane	50 1,2-Dibromo-3-chloropropane
10 1,1-Dichloropropene	23 Tetrachloroethylene	37 n-Propylbenzene	51 1,2,3-Trichlorobenzene
11 Carbon tetrachloride	24 Dibromochloromethane	38 Bromobenzene	52 Hexachloro-1,3-butadiene (Hexachlorobutadiene)
12 1,2-Dichloroethane	25 1,2-Dibromoethane (EDB)	39 1,3,5-Trimethylbenzene	53 Naphthalene
13 Benzene	26 Chlorobenzene	40 2-Chlorotoluene	54 1,2,3-Trichlorobenzene
	27 1,1,1,2-Tetrachloroethane	41 4-Chlorotoluene	

Comparison Direct Injection vs ITEX | Diesel Oil



1 μ l of a Diesel dissolved in Methanol (500ng / μ l) was injected (splitless 15sec.) into the injector to determine a "100% value". 1 μ l of the same solution was added to 12ml water in a 20ml Headspace vial and then analysed with ITEX.

ITEX Parameter:

Extraction Speed: 120 μ l/sec.

Extraction Strokes: 120

Temperature Pumping Syringe / Sample Incubation: 50°C - 10min.

Desorption at 250°C, 15sec. splitless

Chromatography:


Injection: Splitless 15sec. at 250°C / Carrier gas: 0.2bar Hydrogen

Column: BGB-1 15m x 0.32mm ID, 1.0 μ m film

Temperature Program: 40°C - 1min. - 10°C/min to 200°C

Detection: FID 250°C

Environmental Application | VOC, BTEX in Water ppt Level



Chromatogram shows BTEX Compounds at a concentration of 50ng/l using 20 Extraction strokes

Monitoring of BTEX and VOC Compounds in Water (EPA Method 502.2)

Key Words: VOC, BTEX, EPA Method 502.2, ITEX

BTEX and VOC Compounds according to EPA Method 502.2 are analysed using ITEX sample preparation technique. Total sample preparation time of less than 15 minutes allows a high sample throughput.

Sample Preparation:

10ml water is filled into 20ml Headspace sample vials. 3g Sodium chloride and 1µl of the internal standard VOC (50ppb Fluorobenzene in Ethanol) is added. After sample conditioning at 60°C during 10 minutes 20 strokes of the headspace are pumped through the ITEXtrap with a velocity of 100µl/sec. The resulting sensitivity is sufficient to obtain the requested detection limit for drinking water of 0.05µg/l.

ITEX Conditions:

Sample Conditioning @ 60°C, 10 min.

Extraction Strokes: 20 x 1ml

Desorption @ 230°C with 1.3ml Headspace 20µl/sec.

Chromatography:

Column: Rtx-502.2, 60m x 0.32mm, 1.8µm film

Carrier Gas: Helium 20psi

Temperature Program: 40°C - 2 min.- to 240°C - 2 min. at 10°C / min.


Precolumn: 1m x 0.32mm deactivated with DPTMDS

Injector: Gerstel KAS3 with septa @ 150°C isothermal

GC: Varian 3300

Detector: Varian Saturn 4D GC/MS/MS

Perfume Application | Softener-Comparison Porapak Filter Technique vs. ITEX



Chromatography (Thermo TraceGC):

Column	Stabilwax 30m x 0.25mmID x 0.25 mm
Oven	35°C - 0.5min. 15°C / min. 50°C - 0min. 5°C / min. 220°C / 1min
SSL	splitless with surge, surge pressure 20kPa/0.4min (0.5ml), split flow 100ml /0.3min
Carrier	He, 1ml constant flow with vacuum compensation

MS conditions (Thermo Trace MS system):

Ionisation mode	El+
Source temperature	230°C
Interface	220°C
Mass	20-350 amu

ITEX conditions:

Incubation temp	40°C
Incubation time	30 min
Syringe temp	45°C
Extraction volume	1000 ml per stroke
Extraction strokes	10
Extraction speed	100 ml/sec
Desorption temp	200°C
Desorption speed	100 ml/sec
Flush time	5 min

Courtesy of: Givaudan Research Company,
CH-8600 Dübendorf, Zürich. Switzerland, H. Koch

Specifications



ITEX adsorption step out of a sample vial

2.5ml with 1/4" 28 UNF fitting

Extraction Speed:

Selectable from 10 μ l / sec. up to 1000 μ l / sec.

Extraction Strokes:

Selectable from 1 - 999

Extraction Volume:

Selectable from 250 μ l - 2500 μ l / 1 stroke

Desorption Temperature:

+5 $^{\circ}$ C above ambient - 350 $^{\circ}$ C selectable in 1 $^{\circ}$ C increments

Desorption Time:

0 - 300 seconds

Pumping Syringe and Trap Cleaning:

Inert gas purging, 30sec. - 600sec.

Heated Pumping Syringe:

+5 $^{\circ}$ C above ambient - 150 $^{\circ}$ C selectable in 1 $^{\circ}$ C increments

Incubator Oven:

6 heated vial positions for 2ml / 10ml / 20ml vials

+5 $^{\circ}$ C above ambient - 200 $^{\circ}$ C selectable in 1 $^{\circ}$ C increments

Agitation:

Interval shaking 250rpm-750rpm, selectable in 1rpm increments

Incubation Time:

Up to 999 minutes selectable in 1 second increments

Conclusion

This new approach, ITEX or In-Tube Extraction, an automated sample adsorption and thermal desorption technique shows a great potential for volatile compounds in various matrices.

The adsorption at the trap does continuously reduce the concentration of analytes in the headspace and disturbs the equilibrium in the vial. The result is a pumping effect of the analytes from the sample matrix to the headspace. This effect is similar as used by the SPME technique but the active control over the pumping strokes allows the user to utilize rather a dynamic headspace as known by the Purge&Trap technique.

The sensitivity level to be reached can be as far down as the ppt level.

Using the flush gas and a user selectable trap heating temperature does allow to clean the trap very efficiently. Blank runs, even after several hundred injections, do not show a carry over effect. We are aware that this statement can not be generalized. It depends on the sample and matrices. The technique is still too young to be able to give a general statement