

COLLECT PAL

PROTECT YOUR SAMPLE

**NO DRIP
NO LOSS**

**Fraction
Collection**




**Metabolic Profiling
DMPK
Radiolabeled Compounds
ADMET
Metabolomics
Impurity Identification**

**No Loss Technology
Collection During Well-to-Well
Movement
Accurate Volume Delivery**

**Plate Seal
Temperature Control
and
UPLC Compatible**



LEAP Technologies, Inc.
PO Box 969
Carrboro, North Carolina 27510
(919) 929-8814 sales@leaptec.com



**Dedicated User-Friendly
Software**

**Collect into 2-mL vials
96-Well and
384-Well Plates**


**Fastest Collection Speed
@3sec/well**

**Flow Rates
0.1-1.5 mL/min
Tubing id**

0.005" - 0.007" PEEK

**24-plate capacity
Up to 9216 Fractions
In One Run
Multiple Injection
Coordination**

**Specifications:
HTC PAL based
30"x24"x30" Space
Mobile, Flexible**



**LEAP Technologies, Inc.
PO Box 969
Carrboro, North Carolina 27510
(919) 929-8814 sales@leaptec.com**