

# Automated Syringe Exchange Gives More Flexibility



## LEAP PAL-xt with XCHANGE® by SGE

- Large and small volume injections in one run
- 9 Syringe sizes: 10ul, 25ul, 50ul, 100ul, 250ul, 500ul, 1mL, 2.5mL, 5mL
- Syringe is exchanged, aligned and locked automatically
- Reagent and internal standards additions
- Liquid-Liquid extraction in 2ml and 10ml vials
- Serial dilution
- MEPS option for in-syringe SPE
- Syringes are available for LC or GC
- Fast wash station or standard wash station
- Also supported by on-board cycles
- Automatic exchange of syringes with wear during large sample batches



[www.leaptec.com](http://www.leaptec.com)  
[info@leaptec.com](mailto:info@leaptec.com)  
800-229-8814  
[www.leapwiki.com](http://www.leapwiki.com)