

# COMBI PAL

## Can your GC autosampler do this ?

### The CombiPAL.

Static Headspace and  
Liquid injections combined  
in one single instrument

Syringe only concept  
No sample loops  
No transfer lines  
No switching valves



Interfaces with any  
major GC/GC-MS system

Switch from  
Headspace to Liquid  
injection mode by  
fast and easy  
syringe change



Top mounted,  
saves laboratory  
bench space

Stand alone operation  
or Win 95/98/NT®  
remote control software

Optimized for  
high resolution GC  
large volume injections  
high sample throughput

Injects into different  
injector types without  
hardware modifications

CTC Analytics has dedicated the last  
10 years to the continued development  
and high reliability of advanced sample  
injection technology. To learn more about  
the unique GC injector Combi PAL or any  
of our LC/LCMS sample injection systems  
contact your CTC Analytics distributor.



Austria +43 18139993  
Belgium +32 10450025  
Czech Republic +420 2531286  
France +33 562475382  
Germany +49 61261686 & +49 23396037  
UK +44 1462442424  
Israel +972 39318983  
Italy +39 0102471295  
Japan +81 334013898  
Netherlands +31 765411800 & +31 181416555  
Sweden +46 18367799  
Switzerland +41 17323131 & +41 613316020  
USA LEAP Technologies Inc. 800 229-8814  
All other countries +41 617658100

...CTC ANALYTICS

Where design meets performance

CTC Analytics AG  
Industriestrasse 20  
CH-4222 Zwingen  
Switzerland  
Tel: +41 61 7658100  
Fax: +41 61 7658199  
E-mail: info@ctc.ch

